

Instructions For Use KT016-IFU

Rev. Date: Dec. 13, 2018

Revision: 6

Page 1 of 3

626 Wilshire Blvd, Suite 410 Los Angeles, CA 90017 - Tel. +1 (408) 580-1396 www.genebiosolution.com

Giemsa Stain Kit (May-Grunwald) (For Bone Marrow)

Description: The Giemsa Stain Kit (May-Grunwald) is intended for use in the visualization of cells present in hematopoietic tissues and certain microorganisms. This kit may be used on formalin-fixed, paraffin-embedded sections.

Nuclei:	Blue/Violet
Cytoplasm	Light Blue
Collagen:	Pale Pink
Muscle Fibers:	Pale Pink
Erythrocytes:	Gray, Yellow or Pink
Rickettsia:	Reddish-Purple
Helicobacter Pylori:	Blue
Mast Cells:	Dark Blue with Red Granules

Uses/Limitations: For In-Vitro Diagnostic use only.
 Histological applications.
 Do not use past expiration date.
 Use caution when handling these reagents.

Control Tissue: Blood Film
 Any well fixed tissue.

Availability/Contents:

<u>Kit Contents</u>	<u>Volume</u>	<u>Storage</u>
May-Grunwald Stock Solution	500 ml	15-30°C
Giemsa Stock Solution	500 ml	15-30°C
Phosphate Buffer Solution, pH 6.8	500 ml	15-30°C

Precautions: This product is a single-use, non-sterile, in vitro diagnostic device .
 Keep away from open flame.
 Avoid contact with skin and eyes.
 Harmful if swallowed.
 Follow all Federal, State, and local regulations regarding disposal.
 Use in chemical fume hood whenever possible.

Preparation of Reagents Prior to Beginning:

1. Prepare Working May-Grunwald Solution by mixing 25ml of May-Grunwald Solution with 25ml of Phosphate Buffer Solution, pH 6.8.
2. Prepare Working Giemsa Solution by mixing 2.5ml of Giemsa Stock Solution with 50ml of Phosphate Buffer Solution, pH 6.8.

Instructions For Use KT016-IFU

Rev. Date: Dec. 13, 2018

Revision: 6

Page 2 of 3

626 Wilshire Blvd, Suite 410 Los Angeles, CA 90017 - Tel. +1 (408) 580-1396 www.genebiosolution.com

Procedure (Standard):

1. Deparaffinize sections if necessary and hydrate to distilled water.
2. Place slide in staining tray and flood with Working May-Grunwald Solution for 6 minutes. Note: Agitate slide occasionally to insure proper staining.
3. Flood slide with Phosphate Buffer Solution, pH 6.8 until no stain runs off.
4. Flood slide with Working Giemsa Solution for 13 minutes. Note: Agitate slide occasionally to insure proper staining.
5. Flood slide with Phosphate Buffer Solution, pH 6.8 until no stain runs off.
6. Allow slide to remain in Phosphate Buffer Solution, pH 6.8 for an additional 3 minutes.
7. Dip slide quickly in distilled water and air dry at room temperature.
8. Dip slide in Xylene or Xylene Substitute.
9. Mount in synthetic resin.

Procedure (Mast Cells):

1. Deparaffinize sections if necessary and hydrate to distilled water.
2. Place slide in staining tray and flood with Working May-Grunwald Solution for 6 minutes. Note: Agitate slide occasionally to insure proper staining.
3. Flood slide with Phosphate Buffer Solution, pH 6.8 until no stain runs off.
4. Flood slide with Working Giemsa Solution for 13 minutes. Note: Agitate slide occasionally to insure proper staining.
5. Flood slide with Phosphate Buffer Solution, pH 6.8 until no stain runs off.
6. Differentiate by dipping slide in Acetic Acid Solution (0.25%) until background is desired intensity.
7. Dip slide 20 times in Phosphate Buffer Solution, pH 6.8.
8. Dip slide quickly in distilled water and air dry at room temperature.
9. Dip slide in Xylene or Xylene Substitute.
10. Mount in synthetic resin.



Gene Biosolution

626 Wilshire Blvd, Suite
410 Los Angeles, CA
90017

Tel: +1 (408) 580-1396