

## Lambda Light Chain [N10/2]

Format	Catalog no.	Pack size	Dilution
Concentrated	GB-c3063AC	0.1, 1.0 mL	1 : 100
Prediluted	GB-p3063AA	6.0mL	Ready to use

**PRODUCT DESCRIPTION -**

Possible applications of this antibody include the detection of leukemias, plasmacytomas, and certain non-Hodgkin's lymphomas by recognizing the lambda light chains of human immunoglobulins. One hallmark of these cancers is the overexpression of one particular class of light chains. About 2:1 is the typical human kappa/lambda ratio. A kappa/lambda ratio more than 10:1 is indicative of cancerous growth when light chain limitation is clearly present.

**INTENDED USE -**

For In Vitro Diagnostic Use. Lambda Light Chain [N10/2] is a mouse monoclonal antibody designed for laboratory application in the qualitative detection of lambda immunoglobulin light chain protein via immunohistochemistry (IHC) in formalin-fixed paraffin-embedded (FFPE) human tissues. The clinical interpretation of any staining or its absence must be supplemented by morphological studies utilizing appropriate controls and assessed in conjunction with the patient's clinical history and other diagnostic tests by a skilled pathologist.

**SUMMARY AND EXPLANATION -**

This antibody identifies lambda light chains of human immunoglobulins, perhaps aiding in the diagnosis of leukemias, plasmacytomas, and specific non-Hodgkin's lymphomas. The predominant characteristic of these malignancies is the limited expression of a singular light chain class. The typical human kappa/lambda ratio is roughly 2:1. The observation of distinct light chain constriction with a kappa/lambda ratio over 10:1 is indicative of malignant growth.

**PRINCIPLE OF PROCEDURE -**

Antigen identification in tissues and cells is a multi-phase immunohistochemistry procedure. The first step attaches the primary antibody to its designated epitope. Following the tagging of the antigen with a primary antibody, a one-, two-, or three-step detection protocol may be utilized. The one-step approach will utilize an enzyme-conjugated polymer that attaches to the main antibody. A two-step approach will involve the addition of a secondary antibody to bind to the primary antibody. An

enzyme-conjugated polymer is subsequently introduced to bind to the secondary antibody. The three-step detection protocol will include the addition of a secondary antibody to bind to the main antibody, succeeded by a linker antibody step to enhance binding efficacy. An enzyme-conjugated polymer is subsequently introduced to connect with the linker antibody. The presence of bound antibodies is demonstrated by a colorimetric response.

SOURCE - Monoclonal mouse

SPECIES REACTIVITY - Human; others untested

CLONE - N10/2

ISOTYPE - Immunoglobulin G1

Request for lot-specific immunoglobulin concentration.

EPITOPE/ANTIGEN - Lambda immunoglobulin light chain

CELLULAR LOCALISATION - Cytoplasmic

POSITIVE TISSUE CONTROL - Tonsil or Bone Marrow

KNOWN APPLICATIONS - Immunohistochemistry 30-40 min. At RT. Staining of formalin-fixed tissues requires heating tissue sections in between pH 7.4 - 9.0 for 45 min at 95°C followed by cooling at room temperature for 20 minutes.

SUPPLIED AS - Buffer with a protein carrier and preservative

STORAGE AND STABILITY -

Store at 2°C to 8°C. The product is stable to the expiration date printed on the label, when stored under these conditions. Do not use after expiration date. Diluted reagents should be used promptly; any remaining reagent should be stored at 2°C to 8°C.

Materials required but not provided -

- 1) Positive tissue control - Tonsil or Bone Marrow
- 2) Negative control tissue (internal or external)
- 3) Microscope slides and coverslips
- 4) Staining jars or baths
- 5) Timer
- 6) Xylene or xylene substitute

- 7) Ethanol/reagent alcohol
- 8) Deionized/distilled water
- 9) Heating equipment/enzyme for tissue pretreatment step
- 10) Detection system
- 11) Chromogen
- 12) Wash buffer
- 13) Hematoxylin
- 14) Antibody diluents
- 15) Peroxide block
- 16) Light microscope
- 17) Mounting medium

#### LIMITATIONS-

The optimum antibody dilution and protocols for a specific application can vary. These include, but are not limited to fixation, heat-retrieval method, incubation times, tissue section thickness and detection kit used. Due to the superior sensitivity of these unique reagents, the recommended incubation times and titers listed are not applicable to other detection systems, as results may vary. The data sheet recommendations and protocols are based on exclusive use of Genebio products. Ultimately, it is the responsibility of the investigator to determine optimal conditions.